

# PrimePrep™ PCR Purification Kit

## Introduction

*PrimePrep*™ PCR Purification Kit offer simple, rapid and cost-effective method for purification from PCR/enzyme reaction mixtures.

The purified DNA can be directly used in ligation, sequencing and other downstream application.

## Kit Components

Cat. No. Reagents	K-7000 (50 prep.)	K-7001 (200 prep.)
Spin column	50 ea	50 ea x 4
Buffer PCR-B	30 ml	120 ml (40 ml x 3)
Buffer PW	10 ml	30 ml (15 ml x 2)
Buffer PE	10 ml	20 ml

---

## Before you begin

---

- **Add ethanol to Buffer PW before use.**

→ Add 40 ml (K-7001: 60 ml) of absolute ethanol before use.

---

## Experimental Protocol

---

1. Add 5 volumes of Buffer PCR-B to 1 volume of the sample and mix well by vortexing.  
If the PCR reaction product is 50 ul, add 250 ul of Buffer PCR-B.
2. Centrifuge the tube briefly at room temperature.
3. Transfer the mixture to a Spin column.
4. Centrifuge for 30 sec ~ 1 min at 13,000 rpm. Discard the flow-through and re-inserting the spin column to the collection tube.
5. Add 700 ul Buffer PW and centrifuge for 30 sec. at 13,000 rpm.  
Discard the flow-through and re-inserting the spin column to the collection tube.
6. Centrifuge for an additional 1 ~ 2 min at 13,000 rpm to remove residual wash buffer.  
Residual ethanol of washing buffer may inhibitor subsequent enzymatic reaction.
7. Transfer the spin column to new 1,5 ml microcentrifuge tube.  
The 1.5 ml microcentrifuge tube is not provided.
8. Add 50 ul of Buffer PE or deionized distilled water to the center of the membrane in the column, let stand for 1 min and centrifuge for 1 min at 13,000 rpm.  
For larger fragment(>5kb), use pre-warmed (70°C) Buffer PE for best efficiency.